

PreserveTech Folding Commode

HCPCS: E0163

PreserveTech™
Stays Cleaner • Lasts Longer



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Versatile 3-in-1 design enables this model to be used as a bedside commode, raised toilet seat or a toilet safety frame.
- Folds up for convenient, compact storage and transport.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 15.5" to 21.5".
- Find added comfort and support with plastic armrests.
- Commode bucket is easy to empty with a convenient carry handle and pail cover to reduce splashes and spills. Use commode liners for even quicker clean up and less mess!
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- The durable steel frame supports up to 350 lbs.
- Quickly and easily opens and folds, just press push buttons to assemble and disassemble in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.
- Drive PreserveTech products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.
- *Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

PRODUCTS	Description
11148-1A	PreserveTech Folding Commode, 1/1V
SPECIFICATIONS	Description
HCPCS	E0163
Seat Depth	14.5"
Seat To Floor Height	15.5" - 22.5"
Seat Width	13"
Warranty	Limited Lifetime
Weight Capacity	350 lbs
Width Between Arms	18"