




## 18" ALUMINUM WHEELCHAIR/LIGHTWEIGHT TRANSPORT CHAIR DUO

[Home](#) / [Categories](#) / [Wheelchairs & Accessories](#) / 18" Aluminum Wheelchair/Lightweight Transport Chair Duo



### 18" Aluminum Wheelchair/Lightweight Transport Chair Duo: MHWTCDD18

- Made from aluminum
- Standard with nylon upholstery
- Can be used as standard self-propelled wheelchair or a transport chair all in one
- Comes with two sets of wheel locks, one for use with wheelchair and one for the use with transport chair
- Detachable desk length armrest
- Folding back rest for easy transport
- Swing away detachable foot plate/rests with release/lock
- Comes with anti-tippers
- Nylon handle wheel rim
- **300 lbs weight capacity**



### PRODUCT DETAILS



Diagram of a manual wheelchair with various adjustable features labeled:

- Detachable Desk Length Armrest
- Folding backrest
- Back support lock
- 1" Bear wheel / quick release axle
- Swivel
- Wheel lock
- 7" Caster
- 8" Caster
- Swing Away Foot release/lock
- Swing Away Foot



844604104029



**MOBB Health Care Ltd.**  
116A Viceroy Rd. Concord,  
ON L4K 2M1 Canada  
Toll: 1-866-238-7190  
Tel: 905-669-0348  
Fax: 905-669-2796  
[www.mobbbhcc.com](http://www.mobbbhcc.com)