

## Universal Cup Holder

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- New clamp allows cup holder to swivel to keep contents level
- Easily attaches to wheeled walkers, rollators, transport chairs/wheelchairs with tubing diameter from 7/8" to 1"
- Special clamp insert can be added or removed as necessary to accommodate different size tubing
- Made with durable, lightweight plastic

### PRODUCTS

### Description

STDS1040S

Cup Holder, Swivel,1BX,6/CS

### SPECIFICATIONS

### Description

Dimensions

3" (W) x 3" (L).