

## Elevate Upright Walker

HCPCS:



The Drive Elevate Upright Walker brings users increased balance and stability while promoting better posture and providing a more natural line of sight than traditional mobility aids. Height-adjustable upholstered forearm pads keep arms supported and level for comfortable use, while large, ten-inch wheels make it perfect for indoor and outdoor use. With a one-handed, side-to-side fold, the Elevate Upright Walker is the stand-up mobility aid thats easy to take on the go.

- Padded armrests and upright handle adjust to user height preferences
- Durable and comfortable padded seat
- Simple one-handed, side-to-side folding
- Includes removable zipper storage bag with strap
- Dual front fork for optimal stability and maneuverability
- Loop lock brakes for added safety and control
- Adjustable length dual handbrake system
- Large 10 multi-terrain front wheels for indoor and outdoor use
- Height adjustable backrest
- 
- 
- 
- 
- 
- 
- 

PRODUCTS	Description
UPWKRBL10	UPRIGHT WALKER, ROLLATOR, 1EA, BLUE
UPWKRGY10	UPRIGHT WALKER, ROLLATOR, 1EA, GREY
SPECIFICATIONS	Description
Dimensions d'expédition du carton	31" (L) x 12" (W) x 16" (H)
Roulettes	Front: 10 / Rear: 8
Plage de hauteur de la poignée	39 - 44
Matériel	Aluminum
dimensions du produit	31.5 (L) x 25.8 (W) x 48 53 (H)
Poids du produit	23lbs.
Garantie	Limited Lifetime
Capacité de poids	300 Lbs
Handle Extension Range	4" - 6"