



## Transit MOBILITY SCOOTER

Electrically foldable scooter with a remote control. Weight capacity of 115kg (250lb). Up to 11 km of battery life. LED front and back lights.

### Available colors



## SPECIFICATION

Model	S21
Weight Limit	115kg (18 stone)
Overall Dimensions (L x W x H)	86 x 42 x 84cm 34" x 16½" x 33"
Folded Dimensions (L x W x H)	73 x 44 x 42cm 28¾" x 17½" x 16½"
Seat Dimensions (L x W x H)	32 x 40 x 38cm 12½" x 16" x 15"
Footplate Length	37cm (14½")
Wheelbase Length	63cm (25")
Turning Radius	78cm (31")
Ground Clearance	3.5cm 1.4"
Drive Wheel	200mm x 50mm (8" x 2") PU Tyre
Front Wheel	180mm x 40mm (7" x 1.6") PU Tyre
Rear Anti-Tip Castor	Standard
Maximum Speed	3.5 mph (6 kph)
Maximum Range	6 miles 10 km
Stability	3° (Dynamic) 6° (Static)
Weight	19 kg (with battery) 17.5kg (without battery)
Charger Type	2A 110/240V Li-Ion Charger
Controller	50A Dynamic
Motor	180W 24V
Battery	Lithium-Ion Battery 24V 8.7 AH (11.5 AH)