

## essential COLLECTION

## VALUE AND COMFORT

The Essential Collection from Pride<sup>®</sup> Power Lift Recliners provides exceptional value and outstanding comfort. Several fabric options will match any décor while delivering rest and worry-free operation.



Easy-to-use  
hand control with  
USB charging port

**ENERGY  
EFFICIENT**  
WITH BUILT-IN  
BATTERY BACKUP  
(BATTERIES REQUIRED)





## WORRY-FREE OPERATION

### FEATURES & BENEFITS

- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Back has access for adjusting the fill for comfort and fit
- Integrated battery backup ensures chair is operational when the power goes out (Batteries not included)
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Lifetime warranty on steel lift frame
- Footrest extension standard (LC-107 only)

### OPTIONS\*

- Cleaning kit
- Non-skid leg levelers

\* Additional charges apply for all options. Extra lead time may apply for some options.

# Essential LC107





MODEL	LC-107
POSITION	Infinite Position Trendelenburg Chaise Lounger 5'4" - 6'0"
SIZE	Medium
WEIGHT CAPACITY	325 lbs.
SEAT-TO-FLOOR	20"
SEAT DEPTH	20"
SEAT WIDTH	20"
TOP OF BACK TO SEAT	27.5"
FULLY PADDED CHAISE	Yes
SEATING TYPE	Sinuous Spring
DISTANCE FROM WALL	26"
OVERALL WIDTH	35"
CHAIR WEIGHT	116 lbs.
WARRANTY	1-Year Limited

CLOUD 9 FABRICS

Beamsville, ON / 888-570-1113 / [pridemobility.ca](http://pridemobility.ca)

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.