



Base brune – Brown base

Table de lit ProBasics, non inclinable

- ✚ Facilement assemblé avec un minimum de matériel
- ✚ Les roulettes pivotantes à roulement à billes facilitent le déplacement de la table
- ✚ Le dessus en bois laminé durable de 30 po de largeur x 15 po de profondeur offre aux patients suffisamment d'espace pour faire ce qu'ils souhaitent.
- ✚ Base stable en « H » avec revêtement en poudre marron
- ✚ Réglage de la hauteur de 27" à 42" pour s'adapter à n'importe quel patient
- ✚ Supporte jusqu'à 25 lb dans la position la plus basse
- ✚ La table de luxe (modèle n ° PB3037) a une limite de poids de 25 lb

ProBasics Overbed Table, non tilt

- ✚ Easily assembled with minimal hardware
- ✚ Swivel ball bearing casters make table easy to move around
- ✚ 30" W x 15" D durable wood grain laminated top gives patients plenty of room to do what they wish
- ✚ Stable "H" base with brown powder coat finish
- ✚ Height adjustment from 27" to 42" to accommodate any patient
- ✚ Supports up to 25 lbs in lowest position
- ✚ The deluxe overbed table (model no. PB3037) has a weight limitation of 25 lbs