




Original Snap Series modular system


The Original Snap mattress is our modifiable foam surface for the client at high risk of experiencing skin breakdown and may require the addition of ROHO air mattress sections. Firm perimeter borders and snaps hold the interchangeable inserts in place.

Weight capacity: 350 lb

Dimensions :

OSS SYS4 - 35 x 78
OSS SYS4 - 36 x 80
OSS SYS4 - 39 x 80
OSS SYS5 - 35 x 78
OSS SYS5 - 35 x 80
OSS SYS5 - 36 x 80
OSS SYS5 - 39 x 80

Original Snap Series modular system 36" x 80"



Original Snap Series modular system 39" x 80"

