



Molift Mover 205

N14050



Specifications common for all sizes			
Battery	NiMH 26.4 V - 2.6 Ah	Lifting speed 165 lbs (inch/sec)	1.6
Battery charge time (hours)	3	Material	Aluminum, steel
Leg height (inch)	4.5	Maximum user weight (lbs)	451
Length (inch)	50.5	Rated performance	50 lifts at: 165 lbs and 20 inch
Lifting height max (inch)	82.5	Service software	Included
Lifting height min (inch)	14.8	Width chassis, legs open inner (inch)	42
Lifting interval (inch)	56.5	Width chassis, outer (inch)	27.5

	Standard chassis, support arms
Item number	M14050
Weight (lbs)	100.5