



Raised Toilet Seat with Armrests



CODE : 081703990	NEW CODE	Maximum User Weight
Raised Toilet Seat with Armrests	09 155 8873	 150 kg  23½ st



The Homecraft Raised Toilet Seat with Arms allows for increased independence when using the toilet by compensating for impaired movement. This makes it a helpful aid for seniors, people with disabilities, and patients post-hip replacement surgery. Also great as a gift for elderly parents or grandparents.

Choose between three possible seat heights for a custom fit, approximately 2.5”, 4”, and 6” (60, 100, and 150mm). The ideal height allows your feet to touch the floor with your hip joint higher than your knee joint. This keeps your hip and knee joints from being bent past 90 degrees, which is important when recovering from replacement surgery.

Molded plastic base has a removable lid, secure lift-up ABS armrests, and sturdy steel base. The seat can be wiped clean with a disinfectant and the lid removes for easy cleaning. Maximum weight capacity is 300 pounds (136 kg).