

Tröst and Tröst Plus Bed



permobil

Specifications

Weight capacity: The safe working load for the Tröst Bed and Tröst Plus Bed is 600 lbs (272 kg) The safe working load is the total weight, including the person, mattress, head and foot boards, side rails, and other accessories evenly distributed on the bed deck.

TRÖST dimensions : 36" X 80"

Weight: The Tröst Bed weighs approximately 265 lb (120 kg)

The Tröst Plus Bed weighs approximately 281 lb (127 kg)

These weights can vary depending on the options and accessories included

Deck Height Range: 7 – 30" (18 – 76 cm)

Maximum Back Angle: 65 degrees

Maximum Hip to Knee Angle: 30 degrees

Electrical: 120V, 60 Hz