


Système de bain PreserveTech™ Aquachair avec bidet

PreserveTech™ Aquachair Bathing System with Bidet

Article #RTL12A004-WH



CARACTÉRISTIQUES ET AVANTAGES

- Redonnez confiance au bain
 - Chaise de douche unique avec bidet et jet de douche permettant un nettoyage complet du corps sans se tenir debout
 - Choisissez les modes de nettoyage avant ou arrière (le nettoyage avant est plus doux que le jet arrière)
 - Accédez facilement au jet de douche depuis le support de l'accoudoir. Le siège ouvert facilite l'accès aux zones difficiles d'accès.
 - Le siège et le dossier rembourrés antidérapants offrent soutien et confort
 - Ajustez la hauteur en appuyant sur un bouton pour obtenir la position assise idéale. Les réglages de la hauteur du siège sur ce modèle vont de X" à X"
 - Les pieds comprennent une pointe antidérapante et antidérapante pour plus de sécurité sur les surfaces mouillées
 - Le cadre en aluminium léger et durable résistant à la rouille supporte jusqu'à 350 lb.
 - Le dossier, le siège et les embouts traités antimicrobiens aident à prévenir les taches et les odeurs dues à l'accumulation de bactéries, de moisissures et de mildiou.
 - Se connecte facilement à l'alimentation en eau : pas de connexion électrique ni de piles nécessaires
-

FEATURES AND BENEFITS

- Restore confidence in bathing.
- Unique shower chair with bidet and shower jet for complete body cleansing without standing.
- Choose front or rear cleaning modes (front cleaning is gentler than the rear jet).
- Easy access to the shower jet from the armrest support. The open seat provides easy access to hard-to-reach areas.
- Non-slip padded seat and backrest provide support and comfort.
- Adjust the height at the touch of a button for the ideal sitting position.
- Seat height adjustments on this model range from X" to X".
- Feet include a non-slip, non-slip tip for added safety on wet surfaces.
- Lightweight, durable, rust-resistant aluminum frame supports up to 350 lb.
- Antimicrobial-treated backrest, seat and end caps help prevent stains and odors from bacteria, mold and mildew buildup.
- Connects easily to the water supply: no electrical connection or batteries required.

SPÉCIFICATION

Garantie	Limité 1 an
Poids du produit	10,89 livres
Matériel	PP + Aluminium
Capacité de poids	350 livres
Couleur	Blanc
Dimensions assemblé	32,28" (L) X 18,11" (L) X 34,65" (H)
Dimensions de la boîte au détail	12,2" (L) X 17,52" (L) X 21,65" (L)
Poids de la boîte au détail	2,43 livres
Largeur entre les bras	17,8"

SPECIFICATION

Guarantee	Limited 1 year
Product weight	10.89 pounds
Material	PP + Aluminum
Weight capacity	350 pounds
Color	White
Assembled dimensions	32.28" (L) X 18.11" (W) X 34.65" (H)
Retail Box Dimensions	12.2" (L) X 17.52" (W) X 21.65" (L)
Retail box weight	2.43 pounds
Width between arms	17.8"