

PreserveTech Raised Toilet Seat with Bidet

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Control panel activates and adjusts water pressure and cleansing modes
- Bidet provides front or female and back cleansing depending on needs
- Retractable bidet nozzle and guard for a more sanitary experience
- Contoured for safe and comfortable elevation
- Antimicrobial* treated raised toilet seat helps minimize the growth of bacterial odors, mold, and fungi on treated surfaces.
- *Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs, or other disease organisms.
- RTL12C005-WH - Ambient water temperature only
- RTL12C006-WH - Choose between warm and ambient water temperatures
- RTL12C005-WH - Easily connects to the water source behind the toilet
- RTL12C006-WH - Easily connects to the water source behind the toilet and under the sink
- RTL12C006-WH - Comes with an inlet cap to plug warm water connection for ambient water use only to accommodate bathroom layouts

PRODUCTS	Description
RTL12C005-WH	RTS w. Bidet, Ambient Water, 1/CV
RTL12C006-WH	RTS w. Bidet, Warm Water, 1/CV
SPECIFICATIONS	Description
Warranty	One Year Limited
Material	HDPE
Weight Capacity	400 lb
Color	White
HCPCS	E0244
CA Prop 65	no