

Configura® Advance

A highly configurable specialist chair for higher dependency users

The Advance meets many needs in its standard configuration, and can be easily configured to accommodate more complex postural requirements by upgrading the backrest. Additional integrated size adjustments in every dimension and an optional electric upgrade make this a highly configurable and easily recyclable chair.



Overall chair width	700 mm (28")
Overall chair depth	1075 mm (42")
Overall chair height	1280 mm (50")
Overall chair depth with legs elevated	1550 mm (61")
Internal seat width	400 mm - 550 mm (16" - 22")
Seat height (from footplate if applicable)	320 mm - 515 mm (13" - 20")
Floor to seat height	620 mm (24")
Internal seat depth	540 mm (21")
Backrest height from the top of seat cushion	660 mm (26")
Backrest angle	90° - 130°
TIS Angle	Manual version: 7° - 37° Electric version: -10° - 37°
Standard seat angle in neutral	7°
Forward tilt angle	10°
Armrest height from seat	200 mm (8")
The leg rest angle	75° - 60°
Footrest angle	90°, 105°, 120°
The seat height from footplate	320 mm - 515 mm (13" - 20")
Clearance from footplate to the floor	40 mm (2") *Not in tilt forward mode
Seat depth	420 - 550 mm (17" - 22")
Seat width	400 - 550 mm (16" - 22")
Footplate depth	300 mm (12")
The distance between the floor and the bottom part of the frame	130 mm (5")