

# ACCESSORIES

HARVY® Brand, “Nev-a-Slip” Tips • Steel Reinforced Inserts • Contains Latex

## CANE TIPS



## SMOOTH STYLE CANE TIPS



## WALKER AND COMMODOE TIPS\*



## QUAD CANE TIPS



| SIZE   | DESCRIPTION                               | INSIDE   | OVERALL               | OVERALL                   | COLOR AVAILABILITY |                |                 |                 | RETAIL |
|--------|-------------------------------------------|----------|-----------------------|---------------------------|--------------------|----------------|-----------------|-----------------|--------|
|        |                                           | DIAMETER | HEIGHT<br>(in inches) | BASE WIDTH<br>(in inches) | TAN<br>ITEM #      | GREY<br>ITEM # | BLACK<br>ITEM # | BROWN<br>ITEM # |        |
| 8-0396 | Cane Tip                                  | 1/2"     | 1-5/8"                | 1-1/8"                    | 8-0396-01          | 8-0396-33      | 8-0396-08       | —               | x      |
| 8-0397 | Cane Tip                                  | 5/8"     | 1-5/8"                | 1-1/8"                    | 8-0397-01          | 8-0397-33      | 8-0397-08       | —               | x      |
| 8-0398 | Cane Tip                                  | 3/4"     | 1-5/8"                | 1-1/4"                    | 8-0398-01          | 8-0398-33      | 8-0398-08       | —               | x      |
| 8-0399 | Cane Tip                                  | 7/8"     | 1-5/8"                | 1-1/2"                    | 8-0399-01          | 8-0399-33      | 8-0399-08       | —               | x      |
| 8-0392 | Cane Tip                                  | 1"       | 1-7/8"                | 1-1/2"                    | 8-0392-01          | 8-0392-33      | 8-0392-08       | —               | x      |
| 8-2700 | Smooth Style Cane Tip                     | 5/8"     | 1-7/8"                | 1-1/2"                    | —                  | 8-2700-33      | 8-2700-08       | 8-2700-44       |        |
| 8-2800 | Smooth Style Cane Tip                     | 3/4"     | 1-7/8"                | 1-1/2"                    | —                  | 8-2800-33      | 8-2800-08       | 8-2800-44       |        |
| 8-4600 | Quad Cane, Small Base Tip                 | 1/2"     | 1-1/8"                | 1-1/4"                    | —                  | 8-4603-33      | 8-4608-08       | —               |        |
| 8-4600 | Quad Cane, Small Base Tip, retail pkg.*** | 1/2"     | 1-1/8"                | 1-1/4"                    | —                  | 8-4603-23      | 8-4608-28       | —               | x      |
| 8-4700 | Quad Cane, Large Base Tip                 | 5/8"     | 1-1/8"                | 1-1/4"                    | —                  | 8-4703-33      | 8-4708-08       | —               |        |
| 8-4700 | Quad Cane, Large Base Tip, retail pkg.*** | 5/8"     | 1-1/8"                | 1-1/4"                    | —                  | 8-4703-23      | 8-4708-28       | —               | x      |
| 8-1600 | Walker/Commode Tip, round bottom          | 1/2"     | 13/16"                | 7/8"                      | —                  | 8-1600-33      | 8-1600-08       | —               |        |
| 8-1700 | Walker/Commode Tip, round bottom          | 5/8"     | 1-1/4"                | 1"                        | —                  | 8-1700-33      | 8-1700-08       | —               |        |
| 8-1800 | Walker/Commode Tip, round bottom          | 3/4"     | 1-1/2"                | 15/16"                    | —                  | 8-1800-33      | 8-1800-08       | —               |        |
| 8-1900 | Walker/Commode Tip, round bottom          | 7/8"     | 1-1/2"                | 15/16"                    | —                  | 8-1900-33      | 8-1900-08       | —               |        |
| 8-2000 | Walker/Commode Tip, round bottom          | 1"       | 1-5/8"                | 13/4"                     | —                  | 8-2000-33      | 8-2000-08       | —               |        |
| 8-2100 | Walker/Commode Tip, round bottom          | 1-1/8"   | 1-3/4"                | 1-3/4"                    | —                  | 8-2100-33      | 8-2100-08       | —               |        |

\*These tips have been specifically designed to withstand severe pressure incurred by the tendency of walker and commode tubing to cut through the rubber base.

\*\*Ask your customer service representative about product codes for 4/set retail packaged items.

\*\*\*4/set retail package.