

Deluxe Bariatric Shower Chair with Cross-Frame Brace

HCPCS: E0245



- Frame has width/depth cross braces attached with aircraft-type rivets
- Durable, blow-molded plastic bench
- Aluminum frame is lightweight, durable and corrosion proof Drainage holes reduce slipping
- Euro-style clip for leg height adjustment Height adjusts in 1" increments (Figure A)
- Weight Capacity 500 lbs
- Made using recycled materials

PRODUCTS

Description

12021KD-1	Bath BenchHD,w/Bk,1/ea
12022KD-1	1.142Bench,HD,w/o.Bk,1/ea

SPECIFICATIONS

Description

Carton Shipping Weight	32.5 lbs. 8.5 lbs. 27 lbs. 8.0 lbs.
Outside Legs	17.5" (W) x 17.25" (D)
Seat Dimensions	20" (W) x 12" (D)
Seat Height	14" - 19"
Warranty	Limited Lifetime
Weight	7.5 lbs. 6.0 lbs.
Weight Capacity	500 lbs

ACCESSORIES

Description

12021C	Horseshoe Clip,Ht Adjust, 12021, 1/ea, EA
12021T	RubberTip,12021,1/ea, EA