



Mermaid (Electric) with Ranger seat

# Mermaid

The Mermaid hoist is a fixed bath-side hoist that enables patients to safely transfer in and out of the bath. The Mermaid is available in both electric and manual versions and can be fixed either to the side or end of the bath. The Mermaid can also be used with either a fixed seat option or the Ranger transporter system.

## Mermaid features & options

- Electric or manual options
- Wooden and concrete floor fixings are supplied as standard
- Optional side or end fit configuration, safety arms, commode seat and electric retro-fit
- Fixed seat or Ranger transporter system
- Emergency stop/emergency raise (electric option)
- Detachable battery pack (electric option)
- Hand control (electric option)



Mermaid (Manual) with fixed seat option

