



LEVIER POUR LE BAIN À BASE MERMAID RANGER

Safe Working Load: 125 kg / 19.6 st	Maximum Height to Base of Seat: 740 mm
Minimum Height to Base of Seat: 140 mm	Maximum Seat Travel: 485 mm
Maximum Seat Width: 450 mm	Maximum Width Between Arms: 395 mm
Maximum Seat Depth: 410 mm	Centre of Seat to Centre of Mast: 430 mm
Maximum Seat Turning Radius: 850 mm	Maximum Clearance from Mast (to allow for seat rotation): 1000 mm
Total Weight - Mermaid Manual: 29 kg	Total Weight - Mermaid Electric: 37 kg
Ranger Transporter Seat Specifications	Maximum Width: 610 mm
Maximum Depth: 560 mm	Maximum Height: 610 mm
Seat Height: 510 mm	Seat: 510 mm