















## Deluxe Bariatric Shower Chair with Cross-Frame Brace

Item #	<b>Product Description</b>	UOM	
12021-4	With Back, Assembled	4/cs	
12021KD-1	With Back, Knock Down	1/cs	
12022-4	Without Back, Assembled	4/cs	
12022KD-1 Without Back, Knock Down 1/cs		1/cs	

- 500 lb. weight capacity.
- Frame has width/depth cross braces attached with aircrafttype rivets.
- Durable, blow-molded plastic bench.
- Aluminum frame is lightweight, durable and corrosion proof. Drainage holes reduce slipping.
- Euro-style clip for leg height adjustment. Height adjusts in 1" increments. (Figure A)

specifications	
Outside Legs	17.5" (W) x 17.25" (D)
Seat Dimensions	20" (W) x 12" (D)
Seat Height	14" - 19"
Weight (12021   12022)	7.5 lbs.   6.0 lbs.
Carton Shipping Weight (12021-4   12021KD-1)	32.5 lbs.   8.5 lbs.
Carton Shipping Weight (12022-4   12022KD-1)	27 lbs.   8.0 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime

## **Deluxe Bariatric Shower Chair with Cross-Frame Brace**

With Tool-free Removable Back

Item #	<b>Product Description</b>	UOM	
12029-1	Assembled	1/cs	
12029-2	Assembled	2/cs	

- 600 lb. weight capacity.
- Frame has width/depth cross braces attached with aircraft-type rivets.
- Large seat with 16.5" depth.
- Seat height adjusts in 1" increments.
- Durable, blow-molded plastic bench.
- Tool-free removable back.
- Aluminum frame is lightweight, durable and corrosion proof.
- Drainage holes reduce slipping.

specifications				
Outside Legs	21.25" (W) x 18.5" (D)			
Seat Dimensions	16.5" (W) x 16.5" (D)			
Seat Height	16.5" - 21.5"			
Weight	9 lbs.			
Carton Shipping Weight (-1   -2)	10   20 lbs.			
Weight Limit	600 lbs.			
Warranty	Limited Lifetime			