



The Invacare 9000 Topaz Wheelchair features a heavy-duty, durable carbon steel frame and comfortable front riggings and arms for proper positioning. For easy mobility and proven stability, this wheelchair boasts gusseted sides, heavy-duty crossbraces and caster journals.

Low-maintenance, durable carbon steel frame

Simple seat-to-floor height adjustments

Heavy-duty package option features 1000 lb. weight capacity

Durable, triple inner liner keeps seat and back from stretching

Heavy-duty arms designed to accommodate weight capacity

**Overall Width:** 28.5", 30.5", 32.5", 34.5", 36.5", 38.5"

**Overall Height:** 35.5" minimum

**Overall Depth:** 30.5" without riggings

**Seat Depth:** 18", 20"

**Seat To Floor Height:** 15.5"-19.5"; 2" increments

**Back Height:** 16"-18"

**Dimensions:** Folded Width: 15"  
Wheel Height: 24"  
Carton Dimensions: 34.5" x 14.75" x 45"

**Product Weight Capacity:** 700 lb., 1000 lb.

**Shipping Product Weight:** 117 lb.

**Warranty:** Limited Warranty - 3 Year