



**ROTEC**  
INTERNATIONAL

[www.rotecbeds.com](http://www.rotecbeds.com)

VARITECH

# VARITECH



## SPECIFICATIONS

- Load capacity: 499 kg (1100 lb)
- Bed weight: 140.5 kg (309 lb)
- Minimum height: 33.3 cm (13 1/8 in.) - Maximum height: 76.8 cm (30 1/4 in.)
- Maximum back angle: 68 degrees - Maximum foot angle: 34 degrees
- Base dimensions: 88.9 cm (35 in.) 134.6 cm (53 in.) x 198.1 cm (78 in.)
- Overall dimensions with half-section cts: 102.5 cm (41 in.) 149.9 cm (59 in.) x 229.6 cm (90 3/8 in.)
- Recommended mattress height: 12.7 cm (5 in.)
- Trapze support included
- Trendelenburg and reverse Trendelenburg positions included

## DETAILS

- Ultra-quiet DC linear motor
- IP54-approved electrical components protected against dust and water ingress
- 6 receptacles for floor-holding rod - Battery for temporary autonomy
- Castors with central synchronized braking system (2 x 2 castors)
- Bed can be moved to any height, even the lowest
- Electrostatic powder-coated metal surfaces
- Manually adjustable foot section angle (choice of three positions: flexion, standard and extension)
- Bed can be split into two parts for easy transportation
- Adjustable bed base width
- Head and foot bumpers
- Affordable adjustable-width bariatric bed
- Bed with the highest load capacity in its class
- Genuine Canadian product

## OPTIONS

- Rotating half-length headboards
- Headboard and footboard in molded composite or melamine
- Trapeze bar (follows the movement of the bed when height adjustment)
- Nursing control
- Integrated digital scale
- Varitech SHD sectional mattress
- Transport cart

