

GUARDIAN™ Two-Button Folding Walkers

- Easy-to-use, push-button mechanisms let walkers fold
- Each side folds independently to allow easy movement through narrow spaces
- Side cross braces on standard walkers provide extra durability
- Push-button, height-adjustable legs

Compare
walkers
side-by-side
on page 35



G07755



MDS864104



MDS864104B

Premium Red•Dot® Walker

Designed to offer maximum stability and longer life. Aluminum, chrome plated design with steel "K" cross brace for maximum durability. Walker maintains tight tolerance between frame and foot pieces for rattle free use. Deluxe grips and leg tips.

Standard Folding Walker

Lightweight aluminum frame with steel cross brace provides durability. Walker maintains tight tolerance between frame and foot pieces for rattle free use. Comfortable hand grips and non-marking tips. Fits users 5'5" to 6'4".

Basic Two-Button Walker

This walker is a great value for the price! Additional steel components makes this walker more durable and sturdy. 400 lb weight capacity accommodates a greater patient population.



G30755P

Premium Easy Care™ Walker

Lightweight aluminum frame with steel cross brace provides durability. Walker maintains tight tolerance between frame and foot pieces for rattle free use. Comfortable hand grips and non-marking tips

One-Button Folding Walkers

- One-button folding mechanism offers quick, convenient operation
- Folding mechanism can be operated with one hand



G30760P

Premium One-Button Folding Walker

Lightweight aluminum frame allows for easy storage. Designed with superior quality features including reinforced tube connectors at high stress points and tight tolerance for rattle free legs. Deluxe grips and leg tips.



MDS86617

Standard One-Button Folding Walker

Ideal for short-term walker users. Aluminum frame with rugged cross bar provides extra strength. Comfortable grips and non-marking leg tips.