



Model no. 9610-4

Invacare® Safeguard Commode

Model no. 9610-4 (gray-coated steel, four per carton)

- Can be used as a raised toilet seat or a toilet safety frame
- Equipped with a splash shield
- Snap-on seat



Model no. 9669

Invacare® Drop-Arm Commode

Model no. 9669 (one per carton)

Features

- Can be used over the toilet with optional splash shield
- Accommodates transfers in restricted areas
- Drop-arm lever operates simply: lift to lock or push to release
- Gray powder coated, 1-inch welded steel tubing
- Snap-on seat and lid

Specifications

	Base width overall	Width between arms	Seat surface width	Seat opening width	Base depth overall	Seat surface depth	Seat opening depth	Height overall min.-max.	Seat-to-floor height min.-max.	Seat-to-arm height	Product weight	Weight capacity	Units/ carton	Limited warranty
9610-4	25"	19"	16"	8 ½"	19"	14"	10"	23" - 29"	16" - 22"	7 ½"	12 lb.	300 lb.	4	lifetime
9669	20 ¼"	17 ½"	14"	8 ¼"	20 ½"	15"	9 ½"	27" - 31 ½"	17" - 21"	8"	22 lb.	250 lb.	1	3 years

Invacare® Pan and Holder

Part no. 13245X000 (pan)

Part no. 1027004 (holder)

- 8-quart pan capacity



Specifications

Part no.	Length overall	Length inside	Width overall	Width inside
13245X000	16 ¾"	14"	10"	9"

Invacare® Pail with Lid

Model no. 6317 (12 per carton)

- 12-quart pail capacity



Invacare® Splash Shield

Model no. 6319 (12 per carton)

