





Premium Raised Toilet Seat with Removable Arms

Item #	Product Description	Product Seat Depth	UOM
12402	Fits Standard Toilets	17"	1/cs
12403	Fits Elongated Toilets	19.25"	1/cs

- Designed for individuals who have difficulty sitting down or getting up from the toilet.
- Allows individual to use existing toilet seat and lid.
- Tool-free arms are 19.75" apart and can be removed or added as needed.
- Arm sleeves and receivers are made of metal.
- Heavy-duty molded plastic construction.
- Lightweight and portable. Easily attaches to toilet bowl.
- · No tools required for installation.
- Easy to clean.
- Two models accommodate standard and elongated toilets.

specifications	
Seat Dimensions (12402)	13.75" (W) x 17" (D) x 3.5" (H)
Seat Dimensions (12403)	13.75" (W) x 19.25" (D) x 3.5" (H)
Width Between Arms	19.75"
Weight	3.5 lbs.
Carton Shipping Weight	4.5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Toilet Safety Frame

With Height and Width Adjustable Arms

Item #	Product Description	UOM	
12001-4	Assembled	4/cs	
12001KD-1	Knocked Down	1/cs	

- Anodized aluminum is sturdy and lightweight.
- Plastic armrests provide additional comfort and support.
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions). (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet.
- Bracket secures safety frame in place. (Figure B)

specifications					
Width Between Arms	(12001-4)	18.5" (smallest)	22.75" (largest)		
Width Between Arms	(12001KD-1)	18" (smallest) 25" (largest)			
Adjustable Arm Heig	ht	25.5" (lowest) 29.5" (highest)			
Carton Shipping Wei	ght (12001-4)		18 lbs.		
Carton Shipping Wei	ght (12001KD-1)		5 lbs.		
Weight Limit			300 lbs.		
Warranty			Limited Lifetime		
accessories and parts					
12001BR-1 To	let Safety Frame	e Bracket	1/ea		