

DYNAMICO

INDOOR

By Ormesa



www.westechhealth.com



201 Growth Parkway Angola, IN 46703

Canadian Suggested Retail Price List in USD
Jun-18

Phone: 260-665-2769

Fax: 260-665-3047

Email: iim@mobility-usa.com

www.mobility-usa.com

Quote: _____

P.O. _____

Bill To

Dealer: _____ Contact: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Fax: _____

Email: _____ Mark For: _____

Ship To Same as Bill To _____

Dealer: _____ Contact: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Fax: _____

Email: _____ Mark For: _____

Dynamico Features and Specifications

- **Frame** - non-toxic varnished chromium-plated steel
- **Black Pelvic Support** - adjustable in height and depth
- **Black Chest Support** - adjustable in height, removable
- **Black Harness** - double thickness cotton fabric with velcro straps, adjustable in height and girth with belts and buckles, easy to remove for washing
- **Wheels** - 3.2" castors with double brake
- **Bumpers** - semi-elastic PVC to protect from impacts

Dynamico

Part Number	Description
____ DYN-1	Extra Small Dynamico Yellow Frame
____ DYN-2	Small Dynamico Yellow Frame
____ DYN-3	Medium Dynamico Yellow Frame
____ DYN-4	Large Dynamico Black Frame
____ DYN-5	Extra Large Dynamico Black Frame



Positioning Options

Part Number	Description
____ DYN809-1	Arm Supports Xs *
____ DYN809-2	Arm Supports Sm & Md *
____ DYN809-4	Arm Supports Lg & Xl *
____ 943	Arm Straps Xs & Sm
____ 943-3	Arm Straps Md, Lg, & Xl
____ DYN813-1	Adjustable Handlebar Xs *
____ DYN813-2	Adjustable Handlebar Sm & Md *
____ DYN813-4	Adjustable Handlebar Lg & Xl *



Client:

PO:

Dynamico Indoor CAD

Positioning Options Continued

Part Number	Description
_____ DYN815-2	Abduction Saddle Sm
_____ DYN815-3	Abduction Saddle Md
_____ DYN815-4	Abduction Saddle Lg



_____ DYN814-1	Leg Divider Xs
_____ DYN814-2	Leg Divider Sm
_____ DYN814-3	Leg Divider Md
_____ DYN814-4	Leg Divider Lg
_____ DYN814-5	Leg Divider Xl

**Accessories**

Part Number	Description
_____ DYN824-1	Transparent Tray Xs, Sm, & Md *
_____ DYN824-4	Transparent Tray Lg & Xl *



_____ DYN810-2	Weight Bars Sm 11 lbs total
_____ DYN810-3	Weight Bars Md 13.5 lbs total
_____ DYN810-4	Weight Bars Lg 16 lbs total
_____ DYN810-5	Weight Bars Xl 19 lbs total



_____ DYN811	Front Wheels with Adjustable Brakes
--------------	--




_____ DYN812	Directional Locks
--------------	--------------------------




* The following accessories can not be used in combination 809, 813, or 824

Client Measurements

A Back height _____ F Lower leg length _____
 B Chest depth _____ G Head height _____
 C Chest width _____ H Arm height _____
 D Seat depth _____ Age: _____
 E Seat width _____ Weight: _____



Dimensions

Measurements in Inches

Dynamico Indoor	X Small	Small	Medium	Large	X Large
Max Capacity	121 lbs	154 lbs	187 lbs	220.5 lbs	253.5 lbs
User Height	29.5 - 39	37 - 51	47 - 61	55 - 68.5	62.5 - 76.5
Unit Weight	26.5 lbs	32 lbs	34 lbs	37.5 lbs	43.5 lbs
Turning Width	39	41	42	46	50
A	17 - 23.5	26 - 32	30 - 40	33.5 - 44	38.5 - 47
B	20.5 - 27.5	30 - 37	35 - 45.5	36 - 49.5	42.5 - 54
C	7	8	10	12	12.5
D	10	11	12	14	16
E	25.5	29.5	33.5	43	47
F	31	33	34.5	37.5	41.5
G	26	26	26	28	30

