



TILT-IN-SPACE RODEO

Standard Specifications

MEASUREMENTS:

INCHES

HCPCS: E1232

- * Plus or minus .5 inches due to fabric
- ** Reduced height when footplates are flipped
- Δ Weight of chair is without front rigging
- ^ One Piece Angle Adjustable Footplate
- ^^ Individual Angle Adjustable Footplates

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

	RD 10	RD 12	RD 14	RD 16
Seat Width*	10	12	14	16
Seat Depth*	7 - 12	9 - 14	11 - 16	13 - 18
Seat Back Height	21	22	24	25
Seat to floor	20	21	22	22
Recline Adjustment	90° or 110°	90° or 110°	90° or 110°	90° or 110°
Tilt Adjustment	5° - 45°	5° - 45°	5° - 45°	5° - 45°
Overall Height	43	43	44	45
Overall Length	30.5	35	35	36
Overall Width	21	23	25	27
Folded Length	33.5	34	34	40
Folded Width	21	23	25	27
Folded Height	17.5	17	18	18
Headrest Extension	6	8	8	8
Shoulder Strap Height	13 - 19	15 - 21	16 - 22	18 - 24
Footplate Adjustment Angle ^	±20°	±20°	±20°	±20°
Seat to Footplate^	9 - 14	9 - 14	11 - 16	13 - 18
	5.25 - 11.25**	5.25 - 11.25**	7.25 - 13.25**	9.25 - 15.25**
Seat to Footplate^^	7 - 13	7 - 13	9.5 - 15.5	12 - 18
	4.5 - 9.25**	4.5 - 9.25**	4.25 - 11**	7.25 - 13.25**
Weight Capacity Regular/Transit	75/66 lbs	75/66 lbs	100/100 lbs	170/140 lbs
Weight of Chair Δ	29 lbs	33 lbs	37 lbs	39 lbs
Push Handle Adjustment	32 - 48.5	34.5 - 51	35 - 52	36 - 53.5



TILT-IN-SPACE RODEO

Standard Specifications

MEASUREMENTS:

CENTIMETERS

HCPCS: E1232

- * Plus or minus .5 inches due to fabric
- ** Reduced height when footplates are flipped
- Δ Weight of chair is without front rigging
- ^ One Piece Angle Adjustable Footplate
- ^^ Individual Angle Adjustable Footplates

Note: Convaid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

	RD 10	RD 12	RD 14	RD 16
Seat Width*	25.5	30.5	35.5	40.5
Seat Depth*	18 - 30.5	23 - 36	28 - 41	33 - 46
Seat Back Height	53	56	61	63.5
Seat to floor	51	53	56	56
Recline Adjustment	90° or 110°	90° or 110°	90° or 110°	90° or 110°
Tilt Adjustment	5° - 45°	5° - 45°	5° - 45°	5° - 45°
Overall Height	109	109	112	114
Overall Length	77.5	89	89	91.5
Overall Width	53	58.5	63.5	68.5
Folded Length	85	86	86	101.5
Folded Width	53	58.5	63.5	68.5
Folded Height	44.5	43	46	46
Headrest Extension	15	20	20	20
Shoulder Strap Height	33 - 48	38 - 53	41 - 56	46 - 61
Footplate Adjustment Angle ^	±20°	±20°	±20°	±20°
Seat to Footplate^	23 - 36	23 - 36	28 - 41	30.5 - 46
	13 - 28.5**	13 - 28.5**	18 - 34**	23.5 - 39**
Seat to Footplate^^	18 - 33	18 - 33	24 - 39.5	30.5 - 46
	11.5 - 23.5**	11.5 - 23.5**	11 - 28**	18.5 - 34**
Weight Capacity Regular/Transit	34/30 kg	34/30 kg	45.5/45.5 kg	77/63.6 kg
Weight of Chair Δ	13 kg	15 kg	17 kg	18 kg
Push Handle Adjustment	81 - 123	88 - 129	89 - 132	93 - 136



BR0034-SPEC-E1232-ENG-02152017 Printed in U.S.A.