



max 400 mm (15½")



max 75 kg (165 lb)








Manatee User Guide

EN 12182
ISO 17966



2018-08 rev.05

support.R82.org 

**MANATEE 0****MANATEE1****MANATEE 2****MANATEE 3**

mm (inch)

	 0	 1	 2	 3
 Seat width	330 (12¾")	330 (12¾")	330 (12¾")	400 (15½")
 Seat depth	220-250 (8½-9¾")	300-370 (11¾-14½")	390-460 (15¼-18")	400-470 (15½-18¾")
 Back height	430 (16¾")	430 (16¾")	580 (22½")	640 (25")
 Lower leg length	280 (11")	280 (11")	350 (13¾")	350 (13¾")
 Seat height	140-430 (5-16¾")	140-430 (5-16¾")	140-430 (5-16¾")	140-430 (5-16¾")
 Total width	420 (16½")	420 (16½")	420 (16½")	490 (19")
 Length, folded	1030-1100 (40-43")	1100-1170 (43-45¾")	1360-1430 (53-55¾")	1440-1510 (56-58¾")
 Height, folded	130 (5")	130 (5")	130 (5")	130 (5")
 Weight	5½ kg (12 lb)	5½ kg (12 lb)	6 kg (13 lb)	7 kg (15.5 lb)
 Max. load	50 kg (110 lb)	50 kg (110 lb)	75 kg (165 lb)	75 kg (165 lb)

mm (inch)**Shower frame - P/N 882100****Bath tub frame - P/N 882101**

Frame width	590 (22¼")	-	435 (17¼")	-
Frame Length	980 (38½")	-	990 - 1110 (39 - 43¾")	-
Height	-	-	270 (10½")	-
Height, minimum	380 (15")	-	-	-
Height, minimum w/chair	760 (30")	-	-	-
Height middle	440 (17¼")	-	-	-
Height middle, w/chair	840 (33")	-	-	-
Height maximum	480 (19")	-	-	-
Height maximum, w/chair	900 (35½")	-	-	-
Frame height, folded	210 (8¼")	-	-	-
Frame weight	8 kg (17¾ lb)	-	2 kg (4½ lb)	-
Max. load	75 kg (165½ lb)	-	75 kg (165½ lb)	-
Required width of corridor	1500 (59")	-	-	-
Minimum turning diameter	3000 (118")	-	-	-