




480 mm (19")





max 110 kg (242½ lb)



Swan User Guide

EN 12182

2018-04-rev.07

**SWAN 2****SWAN 3****SWAN 4****SWAN 5****SWAN 6****SWAN 7**

mm (inch)	2	3	4	5	6	7
Seat width, back (A)	200 (7¾")	230 (9")	250 (9¾")	300 (11¾")	340 (13½")	410 (16")
Seat width, front (A1)	300 (11¾")	340 (13½")	360 (14¼")	380 (15")	430 (17")	480 (19")
Seat depth (B)	270 (10½")	325 (12¾")	325 (12¾")	380 (15")	460 (18")	460 (18")
Back height (C)	380 (15")	400 (15¾")	420 (16½")	480 (19")	570 (22½")	570 (22½")
Seat height, w/ wheels (E)	370-530 (14½"-20¾")	480-640 (19"-25¼")	480-640 (19"-25¼")	480-640 (19"-25¼")	480-640 (19"-25¼")	480-640 (19"-25¼")
Seat height, w.o./ wheels (E1)	280-440 (11"-17¼")	390-550 (15½"-21¾")	390-550 (15½"-21¾")	390-550 (15½"-21¾")	390-550 (15½"-21¾")	390-550 (15½"-21¾")
Back width, shoulders (F)	225 (8¾")	250 (9¾")	310 (12¼")	390 (15¼")	460 (18")	490 (19¼")
Back width, hips (F1)	240 (9½")	255 (10")	280 (11")	310 (12¼")	355 (14")	430 (17")
Width, standard with wheels (G)	495 (19½")	535 (21")	535 (21")	535 (21")	575 (22¾")	655 (25¾")

mm (inch)	2	3	4	5	6	7
Width, without wheels (G1)	460 (18")	500 (19¾")	500 (19¾")	500 (19¾")	540 (21¼")	620 (24½")
Length, standard with wheels (H)	460 (18")	570 (22½")	570 (22½")	570 (22½")	620 (24½")	620 (24½")
Length, without wheels (H1)	390 (15¼")	500 (19¾")	500 (19¾")	500 (19¾")	550 (21½")	550 (21½")
Height, standard with wheels (I)	680-840 (26¾-33")	710-870 (28-34¼")	780-940 (30¼-37")	790-950 (31-37½")	860-1020 (33¾-40¼")	860-1020 (33¾-40¼")
Height, without wheels (I1)	590-750 (23¼-29½")	620-780 (24½-30¾")	690-850 (27¼-33¾")	700-860 (27½-33¾")	770-930 (30¼-36½")	770-930 (30¼-36½")
Weight	8 kg (17½ lb)	8,5 kg (18¾ lb)	9,5 kg (21 lb)	9,5 kg (21 lb)	10,5 kg (23¼ lb)	11,5 kg (25¼ lb)
Max load	30 kg (66¼ lb)	60 kg (132¼ lb)	70 kg (154¼ lb)	80 kg (176¼ lb)	90 kg (198½ lb)	110 kg (242½ lb)
Aperture	105*155 (4¼*6")	120*195 (4¾*7¾")	125*190 (5*7½")	155*250 (6*9¾")	155*310 (6*12¼")	155*325 (6*12¾")
Angle	18°	18	20°	20°	22°	27°