

# Bariatric Products

BY AMG



**Foot Print (L x W)**  
28" x 20"  
71.1 cm x 50.8 cm

## Bariatric Drop Arm Commode

770-322



This practical institutional commode is designed to facilitate access from any side while providing safety and support to heavier patients. The Extra Wide frame is constructed of a sturdy rust-resistant powder coated steel and the components are constructed of sturdy PVC plastic; for a weight capacity of up to 700 lb. The sturdy arm mechanism easily snaps into place and disengages with a "one finger movement" to make patient transfers easy and safe. It's adaptable, 5 height positions, and its clearance level assures that it can be easily used over a toilet or simply with it's pail. It can also be used with or without the back support to minimize space usage. It features an extra wide width between armrests and extra wide seat; full-length plastic armrests for added support for enhanced patient safety and comfort.

**Pail:** 10 1/2" diameter x 6 1/4" depth

**Easy to use:** Easy, tool-free, pushbutton assembly.

**Superior Compatibility:** Adjustable Leg / Seat Height from 19 3/8" to 23 1/4" (in 1" increments) to accommodate many toilet heights




\*See table for dimensions  
\*Spare parts available



## Technical Specifications

<b>Weight</b>	21.85 lb (9.91 kg)
<b>Clearance</b>	17" (43.2 cm)    21 1/4" (54 cm)
<b>Seat to Floor Height</b>	19 1/4" (48.9 cm)    23 1/2" (59.7 cm)
<b>Overall Height</b>	28" (71.2 cm)    32" (81.2 cm)
<b>Seat to Arm Height</b>	8" (20.3 cm)
<b>Seat Depth</b>	18" (45.8 cm)
<b>Seat Width</b>	23 1/8" (58.7 cm)
<b>Width Between Arms</b>	25 1/4" (64 cm)
<b>Weight Capacity</b>	700 lb (318 kg)

### Dimensions Legend



To learn more about these products or any other AMG product, contact your AMG Representative, call toll-free at 1-800-363-2381 or visit our web site.



AMG Medical Inc.  
Montreal, QC H4T 1V5  
[www.amgmedical.com](http://www.amgmedical.com)

REV. 1.0 1107