



10410

10412



Forearm Crutches, Euro Style, Lightweight Aluminum

Item #	Color	Product Size	Patient Height	Product Handle to Floor Height	UOM
10410	Silver	Adult	5'0" - 6'2"	29" - 38"	1pr/bx
10412	Blue	Adult	5'0" - 6'2"	29" - 38"	1pr/bx

- One-piece molded plastic cuff and hand grip assembly provides safety and comfort.
- Available in attractive colors.

specifications

Weight	1.5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Steel Forearm Crutches

Item #	Product Size	Patient Height	Product Handle to Floor Height	UOM
10403	Adult	5'0" - 6'2"	28" - 37"	1pr/bx
10405	Tall Adult	5'10" - 6'6"	33" - 42"	1pr/bx

- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability.
- Leg and forearm sections adjust independently for optimal sizing.
- Vinyl hand grips are comfortable and durable.
- Extra-large tips for added stability.
- Ortho K Handle is hard plastic.

specifications

Weight Limit	300 lbs.
Weight	2.5 - 2.6 lbs.
Warranty	Limited Lifetime