

DOLPHIN BATH CHAIR SPECIFICATIONS

Model Number	DBC-2081
Width:	31½" (80 cm)
Length (with foot plate):	43" (109 cm)
Overall Height:	Lowered 48.4" (123 cm); Raised 74" (188 cm)
Seat Height (lowered):	20.75" (52.7 cm)
Clearance (seat raised):	43.5" (110.5 cm)
Lift Range:	25.6" (65 cm)
Maximum Capacity:	400 lbs (181 kg)
Front Casters (dual wheel):	3" (8 cm)
Rear Casters (single/locking):	4" (10 cm)
Frame Construction:	Zinc plated stainless steel tubing with powder Coat. Stainless steel chair frame and chair arms.
Chair Construction:	Vinyl coated closed-cell foam with 3/8" PVC plastic core
Outer Mast:	Mounted on durable nylon (plastic) glides for smooth operation
Power System:	<i>See Interchangeable Battery System Technical Data</i>

