

# Molift Smart 150

LIGHT, FOLDABLE AND TRANSPORTABLE WITHOUT TOOLS



## Compact smartness

Molift Smart 150 is designed to be easily wheeled and transported to where it is required. It is designed for home care, for travelling and easy storage. The Molift Smart is foldable in one piece without using any tools.

## Low weight and unique construction

Due to its low weight, the small footprint and the unique construction, Molift Smart 150 is easy to manoeuvre even in narrow spaces.

## Comfortable and safe transfer

The hoisting range of 27–168 cm (10.6"–66.1") provides an excellent maximum hoisting height and at the same time gets very low, facilitating hoisting from the floor. Standard 4-point sling suspension ensures a comfortable and spacious hoisting position for the user.

Molift Smart 150 has been awarded for Design Excellence by the Norwegian Design Council.

|  |                          |  |  |   |
|--|--------------------------|--|--|---|
| <b>Description</b><br>Molift Smart 150   | <b>Item no.</b><br>10200 | <b>Leg height</b><br>110 mm (4.4")                                       | <b>Hoisting interval</b><br>1410 mm (55.5")                                      | <b>Cleaning instructions</b><br>Use pH-neutral cleaning agents. If disinfection is required, use isopropyl alcohol. Do not use abrasive or corrosive substances for cleaning. |
| <b>Service Software</b><br>Included  |                          | <b>Width chassis/legs outer</b><br>660 mm (25.9")                        | <b>Battery</b><br>NiMH 14.4 V – 2.2 Ah<br>Item no. 0541000                       |   |
| <b>Material</b><br>Aluminium, steel  |                          | <b>Width open legs inner</b><br>925 mm (36.4")                           | <b>Battery charge time</b><br>Ca 3 hours   |   |
| <b>Total Weight</b><br>26 kg (57 lbs)  |                          | <b>Dimension folded</b><br>1160 x 475 x 360 mm<br>(45.32" x 18.7" x 14") | <b>Hoisting speed</b><br>60 mm/sec (2.4"/sec)                                    |   |
| <b>Weight when separated</b><br>Part 1: 13 kg (28.6 lbs)<br>Part 2: 13 kg (28.6 lbs) |                          | <b>Hoisting height</b><br>270–1665 mm<br>(10.6"–66.1")                   | <b>Rated performance</b><br>40 hoists at: 75 kg and 50 cm<br>(165 lbs and 19.5") |   |
| <b>Length</b><br>1200 mm (47.2")   |                          |  |  |   |



## Features



### 4-point suspension

Standard 4-point suspension ensures comfortable and more spacious hoisting position for the user.



### Hoists easily from floor

Molift Smart has a large hoisting span that allows hoisting from the floor as well as high surfaces.



### Foldable without tools

Designed for homecare and travelling.

## Accessories



### Travel Suitcase – Molift Smart 150

This hard shell suitcase provides additional protection when transporting by air. It is TSA approved for air travel and ideal for the client with an active life.  
Item no. 0990200



### Travel Bag for Molift Smart 150

The soft shell travel bag has extra large transport wheels. It protects the hoist, and enables easy transporting.  
Item no. 3049175



Molift Smart folds and unfolds easily without using any tools and can be separated into two parts where the heaviest part weighs only 13 kg.