



Model no. 12120-37-746

Dolomite Maxi+ Rollator

Basket and back support standard.

Model no. 12120-37-746      MAXI+ Standard  
Model no. 12122-37-746      MAXI+ Low  
Colour: Navy Blue

Dolomite’s Maxi+ rollator offers high quality and function for users requiring additional seating space and increased weight capacity. The Maxi+ has a reinforced frame with a maximum weight capacity of 440 lbs.

Features

- Forward swept handles to promote better posture
- Flip-up seat for greater reach access and gait training
- Anatomical handles
- Patented braking system
- Extra wide seating area
- Large 8" SLII wheels with superior bearings for increased maneuverability

Specifications

	Low	Standard
Distance between rear wheels	22"	22"
Minimum handle height	29"	32"
Maximum handle height	35"	38 ½"
Seat height	22"	25 ½"
Overall width	27 ¼"	27 ¼"
Overall length	25 ½"	29"
Width between handles	21 ½"	21 ½"
Rollator weight (lbs.)	22	22
Weight capacity (lbs.)	440	440