

# Etac Relax shower seat

*Scandinavian design for the modern bathroom*

## Attention to detail

Etac Relax has been designed to ensure a secure and enjoyable shower experience. Design, function and safety are extremely important to us, and this is presented here in the form of a discreet, wall-mounted shower seat. The shape of the seat offers the option of varied sitting positions, resulting in improved access. Relax is always close at hand, but as it can be folded upwards it never takes up any unnecessary space.

## Adapted to your needs

Etac Relax comes in two colours and six different models, which have all been designed to create options specifically for you. The shower seat can easily be supplemented with accessories like supporting legs or back support.



### Product name

- White	81703010
- White, with arm supports	81703020
- White, with supporting legs	81703030
- White, with arm supports and back support	81703060
- White, with arm supports and supporting legs	81703040
- White, with arm supports, back support and supporting legs	81703050
- Volcano grey	81708000
- Volcano grey with arm supports	81708010
- Volcano grey with supporting legs	81708020
- Volcano grey with arm supports and back support	81708050
- Volcano grey with arm supports and supporting legs	81708030
- Volcano grey with arm supports, back support and supporting legs	81708040

### Colour

White, Volcano grey

### Maximum user weight

With/without supporting legs: 150/125 kg

### Material

Seat; polypropylene

Arm supports; polypropylene and TPE

Legs; aluminium and stainless steel,

Ferrules; TPE



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F if so required.