



Lo-Back Contour Wrap-Around in Stingray Gray  
Available in 3 sizes

#### FEATURES

- Lo-Back™ and Hi-Back™ options available (Hi-Back model includes four-point harness for added trunk support)
- Bacteria and fungal resistant ComfortMesh™ seat
- Latex and silicone free



Hi-Back Contour Wrap-Around in Stingray Gray  
Available in 3 sizes

#### BENEFITS

- Orthopedic back height adjusts vertically
- Wide-stance stable base fits into any tub
- Easily adjustable back heights
- Compact design for easy storage or travel
- Rustproof frame and tear resistant mesh fabric for extra durability
- Mounting Base allows you to upgrade or replace the support back to the next size for a fraction of the cost of a new support

#### ACCESSORIES



Padded Abductor



Pelvic Support Belt

## ORDER FORM FAX TO 905-787-8550 OR SUBMIT ONLINE

Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ E-mail: \_\_\_\_\_

Ship to address: \_\_\_\_\_

BASE UNIT	LO-BACK WRAP-AROUND			HI-BACK WRAP-AROUND		
	BCWL-7210S	BCWL-7310M	BCWL-7410L	BCWH-7610S	BCWH-7710M	BCWH-7910L
Product Size	Small	Medium	Large	Small	Medium	Large
Weight Capacity	N/A	N/A	N/A	N/A	N/A	N/A
Overall Depth	29"	29"	29"	29"	29"	29"
Seat Depth	14.5"	14.5"	14.5"	14.5"	14.5"	14.5"
Overall Seat Width	18"	18"	18"	18"	18"	18"
Inside Seat Width	15"	15"	15"	15"	15"	15"
Inside Back Width at Top	7"	11"	13"	7"	8"	11"
Inside Back Width at Bottom	7"	11"	13"	12"	12"	14"
Back Height Range	9" - 14"	10" - 14"	11" - 16"	17" - 22"	19" - 24"	24" - 28"
Product Weight	5.5 lb	6.5 lb	7.75 lb	9.5 lb	9.75 lb	12 lb

### ACCESSORIES (fit all sizes of bath chair unless noted)

Padded Abductor BCWZ-A7815

Pelvic Positioning Belt BCWZ-A7830