


NIGHTHAWK

SPECIALTY MADE BACKS

ADJUSTABLE BACKS

CANADA'S #1 Adjustable Back

The Adjustable Back uses five strips of webbing and velcro for user friendly adjustments. The top padding covers over strapping both front and back. Nighthawk is able to manufacture these backs to fit most models of wheelchairs. For added privacy, dignity flaps are available at an additional cost.




Velcro Around Cane



Slide Over Cane



Screws On To Cane



Heavy Duty




Standard Foam



Egg Crate Foam With Lycra



Gel



Pin Core Foam

Benefits:

- Those With Fixed Postural Deformities
- Relieve Pressure Caused By Boney Prominences
- Accommodates Clients With Different Heights And Sizes
- Reduce The Amount Of Shearing
- Provide Privacy Between The Back And Seat
- Relieve Discomfort From Sitting For Long Periods Of Time
- Adjustable Support For Different Areas As Fatigue Sets In Throughout The Day

Standard Padding



PART# DESCRIPTION

2026-B There is no need to remove back canes in order to attach this back, simply velcro on to back canes and use top grommet hole to attach upholstery with screws.

DESCRIPTION

2026-C This back must slip down the back canes in order to mount. Grommet holes are available to attach upholstery with screws.

DESCRIPTION

2026-D This back mounts with rails and screws.
Provide Hole Pattern

PART# DESCRIPTION

2026-B/2148 There is no need to remove back canes in order to attach this back, simply velcro on to back canes and use top grommet hole to attach upholstery with screws.

DESCRIPTION

2026-C/2148 This back must slip down the back canes in order to mount. Grommet holes are available to attach upholstery with screws.

DESCRIPTION

2026-D/2148 This back mounts with rails and screws.
Provide Hole Pattern

Egg Crate Foam & Lycra



Spine Gel Back



PART# DESCRIPTION

2026-B/2143 There is no need to remove back canes in order to attach this back, simply velcro on to back canes and use top grommet hole to attach upholstery with screws.

DESCRIPTION

2026-C/2143 This back must slip down the back canes in order to mount. Grommet holes are available to attach upholstery with screws.

DESCRIPTION

2026-D/2143 This back mounts with rails and screws.
Provide Hole Pattern

NIGHTHAWK

SPECIALTY MADE BACKS

ADJUSTABLE BACKS

Pincore Foam

NEW



PART# DESCRIPTION

2026-B/2151 There is no need to remove back canes in order to attach this back, simply velcro on to back canes and use top grommet hole to attach upholstery with screws.

DESCRIPTION

2026-C/2151 This back must slip down the back canes in order to mount. Grommet holes are available to attach upholstery with screws.

DESCRIPTION

2026-D/2151 This back mounts with rails and screws.

Provide Hole Pattern



Heavy Duty Standard Foam



PART# DESCRIPTION

2026HD Instead of The usual 5 straps at 1", We will incorporate 3 straps at 2" with 2 straps at 1" for extra support. Available only in slide over and screw on styles for safety.
From 21" - 24" Wide

Xtra Heavy Duty Standard Foam



PART# DESCRIPTION

2026XHD Instead of The usual 5 straps at 1", We will incorporate 3 straps at 2" with 2 straps at 1" for extra support. Available only in slide over and screw on styles for safety.
From 25" - 30" Wide


NIGHTHAWK

SPECIALTY MADE BACKS

ADJUSTABLE BACKS


Replacement Adjustable Back Pads

Standard Padding




PART#	DESCRIPTION
2026-REPSF	Replacement Adjustable Back Pad Standard Foam

Egg Crate Foam & Lycra




PART#	DESCRIPTION
2148-REPEC	Replacement Adjustable Back Pad Egg Crate Foam

Pincore Foam



PART#	DESCRIPTION
2151-REPPC	Replacement Adjustable Back Pad Pincore Foam

Ovation® Gel
Spine Gel Back



PART#	DESCRIPTION
2143-REPOG	Replacement Adjustable Back Pad Ovation Gel & Standard Foam

Back Cover Only

Call Us For Custom Made Adjustable Back Pads