

WE HAVE A SIX ACRE PROPERTY AND THE EXTREME X8 **ALLOWS ME TO GET** TO ALL THE UNEVEN TERRAINS AND SEE PARTS OF THE FARM THAT I HAVE NOT BEEN ABLE TO ACCESS IN FOUR YEARS! THANK YOU FOR GIVING ME MY FREEDOM AND INDEPENDENCE."

Michael Trost Maryland USA

HIT THE ROAD

Or the sand, or the bush, or the snow.*

*Not to be driven through water

MAKING YOUR CHAIR **WORK FOR YOU**

Magic Mobility wheelchairs can be found on building sites, farms, stables, national parks all allowing people to continue to work.



DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help you find them, challenge yourself and embrace your freedom.





Exclusive USA distributor since 1999

201 Growth Parkway Angola, IN 46703 **United States**

Tel 1-260-665-2769 Email iim@mobility-usa.com Web mobility-usa.com





magic Mobility

Extreme





MAGIC MOBILITY EXTREME X8

The Magic Mobility Extreme X8 gives you the freedom that you always knew you had. The X8 will be your partner to just about any destination. Sand, snow and every terrain you can throw at it. With unrivaled off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

Basic specifications:**

- Length: 451/2"
- Width: 28"
- Seat to floor height: 18½" (19½" with swingaway legrests)
- Weight limit: 400lbs
 340lbs with lift and tilt
- Top speed: 6.2mph
- Motors: 700W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Batteries: Group 24, 75Ah
 LA/AGM and Gel
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-23" wide).

Rehab positioning:

- Power options: 12" seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest and adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power elevating
- Footplates: 1 or 2 piece fixed/ flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, gun mounts, bags, pockets, keep safe secret boxes, tables, drinks holder, oxy and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Crash test compliant tie downs
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:*

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living.

Options:



SKY BLUE



BLACK



TANGERINE



specifications subject to change. Check online for current version

^{**}Various model options excluded





Magic Mobility +61 3 8791 5600 sales@magicmobility.com.au

Explore with us magicmobility.com.au













"WE'VE BEEN

CREATING CUSTOM MADE VEHICLES FOR 25 YEARS, SO


THAT YOU CAN LIVE

THE LIFE YOU CHOOSE"

Ashley Daff

Director

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.