

Two Wheeled Walker with Seat

HCPCS: E0143/E0156-1,E0143+E0156



- Features comfortable and convenient padded seat and backrest
- Quickly folds when not in use
- Rear leg tips act as brakes when pressed down
- Height adjustable handles
- Large swivel wheels for added maneuverability
- Contoured handles for a comfortable grip

PRODUCTS	Description
1239RD	2 Wheeled Walker w/Seat, Red

SPECIFICATIONS	Description
Carton Shipping Weight	14 lbs
Handle Height	29" - 38"
Overall Width	22"
Seat Dimensions	15.875" (W) x 8" (D)
Seat Height	20.5"
Warranty	Limited Lifetime
Weight	12 lbs
Weight Capacity	300 lbs

PARTS	Description
9501S72603	Backrest Pivot Brackets for R726,1 pr, PR
95032F260N10	Height Adjustment Knobs,260N, 10257 1pr, PR
95032F260N06	Dust Cap,260N,10257, 1/pr, PR
1239-09	Handlebar W/Hand Grip, Left, 1239, 1/ea, EA
1239-10	Handlebar W/Hand Grip, Right, 1239, 1/e, EA
1239-07	Backrest, 1239, 1/ea, EA
1239-06	Seat W/Brackets & Screws, 1239, 1/ea, EA
1239-02	Caster, 1239, 1/ea, EA
1239-05	Frame Plug, 1239, 1/ea, EA
1239-03	Folding Hinge W/Hardware, 1239, 1/ea, EA
1239-01	Rear Glide, 1239, 1/ea, EA
1239-08	Backrest Release Pin, 1239, 1/ea, EA
1239-04bl	Fork, 1239BL, 1/ea, EA