

# Don't Let Them Walk Away!

New & Affordable!

## Prevent Wandering At Home With Anti-Wandering DepartAlert®




### All-in-One

This "All in One" Anti-Wandering Door Bar System functions with just a door alert, contact switch, and resident wristband. Simply mount the door bar by any exit or doorway in your home or facility, plug it in, attach contact switch and place wristband on resident. When a resident wearing a wristband opens the door bar or attempts to go through the doorway, the door bar will alert audibly and visually.

### Simple Installation

1. Mount Door Bar
2. Activate Wristband Transmitter
3. Test



### Affordable – No Wiring – No Electricians

MFR#	PRODUCT DESCRIPTION
TL-4005HOME	Single Door Bar Monitor System (1 year warranty)
TL-4012	Resident Wristband Transmitter (6 month warranty)
TL-4014W	Wall Mount Reset Button (1 year warranty)

### Optional Additional Components:



TL-2025  
Magnetic  
Door Lock



TL-4015  
Central Monitor



5 Extra Bands

\*Facility Door monitoring also available with the TL-4005SYS

Anti-Wandering Solutions  
**SMART**  
Quieter Fall Prevention™

Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)