

# Configura® Comfort

## The most versatile riser recliner

Building on the Lite 2, the Comfort is a dual motor tilt-in-space chair with independent no-shear backrest recline, offering standard and flat rise options. The subframe and slide back armrests are also designed to be compatible with a range of transfer aids.



	Extra Small	Small	Small Tall	Medium	Medium Low Profile	Large	Large Low Profile
Overall chair width	700 mm (27.5")	750 mm (29.5")	750 mm (29.5")	800 mm (31.5")	800 mm (31.5")	850 mm (33.5")	850 mm (33.5")
Overall height	1080 mm (42.5")	1080 mm (42.5")	1110 mm (43.7")	1110 mm (43.7")	1080 mm (42.5")	1115 mm (43.9")	1080 mm (42.5")
Overall length	885 mm (34.8")	885 mm (34.8")	930 mm (36")	930 mm (36")	930 mm (36")	940 mm (37")	940 mm (37")
Overall chair depth with legs elevated	1565 mm (61.6")	1565 mm (61.6")	1660 mm (65.4")	1660 mm (65.4")	1660 mm (65.4")	1660 mm (65.4")	1660 mm (65.4")
Internal seat width	400 mm (16")	450 mm (18")	450 mm (18")	500 mm (20")	500 mm (20")	550 mm (22")	550 mm (22")
Floor to seat height	400 mm (16")	400 mm (16")	450 mm (18")	450 mm (18")	400 mm (16")	450 mm (18")	400 mm (16")
Internal seat depth Position 1	450 mm (18")	450 mm (18")	450 mm (18")	450 mm (18")	450 mm (18")	450 mm (18")	450 mm (18")
Internal seat depth Position 2	500 mm (20")	500 mm (20")	500 mm (20")	500 mm (20")	500 mm (20")	500 mm (20")	500 mm (20")
Backrest height from seat	710 mm (28")	710 mm (28")	710 mm (28")	710 mm (28")	710 mm (28")	710 mm (28")	710 mm (28")
Chair to wall/furniture	570 mm (22.4")	570 mm (22.4")	630 mm (24.8")	630 mm (24.8")	630 mm (24.8")	605 mm (23.8")	605 mm (23.8")
Backrest angle – std	100°	100°	100°	100°	100°	100°	100°
Backrest angle – recline	128°	128°	128°	128°	128°	128°	128°
Rise angle – std	23°	23°	23°	23°	23°	23°	23°
Rise angle – med rise	14°	14°	14°	14°	14°	14°	14°
Rise angle – flat	2°	2°	2°	2°	2°	2°	2°
Rise height at front of seat – std	510 mm (20.1")	510 mm (20.1")	545 mm (21.5")	545 mm (21.5")	510 mm (20.1")	545 mm (21.5")	510 mm (20.1")
Rise height at front of seat – med	570 mm (22.4")	570 mm (22.4")	620 mm (24.4")	620 mm (24.4")	570 mm (22.4")	620 mm (24.4")	570 mm (22.4")
Rise height at front of seat – flat	660 mm (26")	660 mm (26")	710 mm (28")	710 mm (28")	660 mm (26")	710 mm (28")	660 mm (26")
TIS angle	26°	26°	26°	26°	26°	26°	26°
Standard seat angle in neutral	6°	6°	6°	6°	6°	6°	6°
Forward tilt angle	n/a	n/a	n/a	n/a	n/a	n/a	n/a
Armrest height	180 mm (7.1")	180 mm (7.1")	180 mm (7.1")	180 mm (7.1")	180 mm (7.1")	180 mm (7.1")	180 mm (7.1")
Weight capacity	160 kg (25 st)	160 kg (25 st)	160 kg (25 st)	160 kg (25 st)	160 kg (25 st)	160 kg (25 st)	160 kg (25 st)
Overall weight	65 kg (143 lbs)	65 kg (143 lbs)	68 kg (150 lbs)	68 kg (150 lbs)	68 kg (150 lbs)	71 kg (156 lbs)	71 kg (156 lbs)