

R82 Crocodile

The R82 Crocodile is a significant advancement in children's walking aids, designed with the potential to support and enhance mobility for children facing diverse challenges. As a child walker for disabled children, it's crafted to potentially foster greater independence and freedom.

The R82 Crocodile is a pediatric walking aid, designed to support and enhance mobility for children with physical limitations.

The Crocodile's adjustable nature and ease of use makes it a reliable gait trainer. Its thoughtful design, safety considerations and adaptability make it an essential assistive device for a wide range of users. The Crocodile can be modified by adding accessories according to the users' individual requirements for postural support. It is available in seven sizes making the Crocodile your option for the child's first steps or stable support later in life. It is the perfect companion for everyday adventures.

The Crocodile sizes 0, 1, 2 and 3 features smaller, slimmer wheels, making it perfect for indoor use or on smooth even outdoor surfaces. These wheels are designed for easy maneuverability, providing smooth and effortless movement. With an adjustable tension on the rear wheels, it offers a reduced rolling resistance and a controllable ride.

The Crocodile AT sizes are better suited for a more rugged environment, with their larger wheels with a bigger footprint and a deeper tread, that excel on surfaces like gravel, dirt paths and in uneven terrains. This ensures a better grip, delivering a smoother ride and better comfort on gravel and grass, while improving stability and traction. This opens more opportunities to challenge postural control and balance and offers additional possibilities for enjoying outdoor playgrounds.



FEATURES

Ready for adventure

The Crocodile is the perfect companion for adventures. It has a lightweight and foldable design making it easy for transportation. Take your Crocodile with you wherever you go!



The right fit

The Crocodile has a height adjustable top frame and multi adjustable handles, along with a wide range of accessories, making it easy to adjust to your child's...


[Read more](#) ↗



All Terrain - all the time

Our Crocodile comes in two size variants; Crocodile and Crocodile All-Terrain with a total of 7 sizes. The AT sizes have larger wheels for an easier manoeuvrability o...

[Read more](#) ↗



Safety in all situations

When exploring the world, safety is a top priority. Therefore, the Crocodile comes with an anti-reverse wheel stop and is available with both hand brakes and fixe...

[Read more](#) ↗

Size 1 AT & 2 AT



Size 3 AT



Size 0, 1 & 2



Size 3





Size 1 AT & 2 AT



Size 3 AT



Size 0, 1 & 2



Size 3



Variants

	Canadian french	Crocodile Taille 1	Size 1 AT	Crocodile Taille 2	Size 2 AT	Crocodile Taille 3	Size 3 AT
[ItemR82ItemNumber]	86800	86801	86801AT	86802	86802AT	86803	86803AT
[ItemBrakeGripDistanceMM]	74	74	74	74	74	74	74
[ItemHandgripLengthMM]	190	190	190	190	190	250	250
[ItemHandgripWidthMM]	25	32	32	32	32	32	32
[ItemHeightFloorToHandlesWstarterkitMM]	-	380-540	-	-	-	-	-
[ItemLargeCastorMM]	W=30/ D=200	W=30/ D=200	W=55/D=290	W=30/ D=200	W=55/D=290	W=30/ D=200	W=55/D=290
[ItemLengthMM]	640	720	780	840	900	870	930
[ItemLengthWstarterkitMM]	-	660	-	-	-	-	-
[ItemMassDeviceKG]	6.1	6.4	8	6.8	8.4	11.3	12.8
[ItemMaxRatedLoadKG]	38.1	38.4	40	53.8	55.4	93.3	94.8
[ItemMaxUserHeightMM]	1050	1150	1190	1500	1540	1800	1840
[ItemMaxUserMassKG]	30	30	30	45	45	80	80
[ItemSmallCastorMM]	W=28/D=140	W=28/D=140	W=50/D=190	W=28/D=140	W=50/D=190	W=28/D=140	W=50/D=190
[ItemTurningDiameterMM]	1440	1600	1600	1900	1900	2000	2000
Folded width (mm)	580	620	660	690	700	690	720
Hauteur du sol aux poignées (mm)	300-610	420-710	460-750	560-850	600-890	740-1025	780-1065
Hauteur, replié (mm)	300	300	345	270	310	290	330
Largeur du châssis (mm)	580	630	660	690	700	690	720
Larghezza tra manopole (mm)	100-250	190-350	190-350	240-420	240-420	330-480	330-480
Longueur, replié (mm)	650	720	750	860	890	1000	1030
Static stability downhill (°)	15	15	15	15	15	15	15
Static stability sideways (°)	4.5	4.5	4.5	4.5	4.5	4.5	4.5
Static stability uphill (°)	7	7	7	7	7	7	7

Size 1 AT, 2 AT & 3 AT				
		mm (inch), kg (lb)		
		Size 1 AT	Size 2 AT	Size 3 AT
	Width between handles	190 - 350 mm (7½" - 13¾")	240 - 420 mm (9½" - 16½")	330 - 480 mm (13" - 19")
	Height from floor to handles	460 - 750 mm (18" - 29½")	600 - 890 mm (23½" - 35")	780 - 1065 mm (30¾" - 42")
	Overall width	660 mm (26")	700 mm (27½")	720 mm (28¼")
	Length	780 mm (30¾")	900 mm (35")	930 mm (36½")
	Turning diameter	1600 mm (63")	1900 mm (75")	2000 mm (78¾")
	Maximum user height	1190 mm (46¾")	1540 mm (60¾")	1840 mm (72½")
	Mass device	8 kg (17.6 lb)	8.4 kg (18.5 lb)	12.8 kg (28.2 lb)
	Maximum user mass	30 kg (66 lb)	45 kg (99 lb)	80 kg (176 lb)

Size 1 AT, 2 AT & 3 AT				
		mm (inch), kg (lb)		
		Size 1 AT	Size 2 AT	Size 3 AT
				
	Static stability uphill	7°	7°	7°
	Static stability sideways	4.5°	4.5°	4.5°
	Folded length	750 mm (29½")	890 mm (35")	1030 mm (40½")
	Folded width	660 mm (26")	700 mm (27½")	720 mm (28¼")
	Folded height	345 mm (13½")	310 mm (12¼")	330 mm (13")

Size 1 AT, 2 AT & 3 AT				
		mm (inch), kg (lb)		
		Size 1 AT 	Size 2 AT 	Size 3 AT 
	Maximum rated load	40 kg (88.1 lb)	55.4 kg (122.1 lb)	94.8 kg (209 lb)
	Small castor (width / diameter)	W=50 mm / D=190 mm (W=2" / D= 7½")	W=50 mm / D=190 mm (W=2" / D= 7½")	W=50 mm / D=190 mm (W=2" / D= 7½")
	Large castor (width / diameter)	W=55 mm / D=290 mm (W=2¼" / D=11½")	W=55 mm / D=290 mm (W=2¼" / D=11½")	W=55 mm / D=290 mm (W=2¼" / D=11½")
	Handgrip width	32 mm (1¼")	32 mm (1¼")	32 mm (1¼")
	Handgrip length	190 mm (7½")	190 mm (7½")	250 mm (9¾")
	Brake grip distance	74 mm (3")	74 mm (3")	74 mm (3")
	Static stability downhill	15°	15°	15°

Size 0, 1, 2 & 3					
		mm (inch), kg (lb)			
		Size 0	Size 1	Size 2	Size 3
	Width between handles	100 - 250 mm (4" - 9¾")	190 - 350 mm (7½" - 13¾")	240 - 420 mm (9½" - 16½")	330 - 480 mm (13" - 19")
	Height from floor to handles	300 - 610 mm (11¾" - 24")	420 - 710 mm (16½" - 28")	560 - 850 mm (22" - 33½")	740 - 1025 mm (29¼" - 40¼")
	Height from floor to handles with starter kit	N/A	380 - 540 mm (15" - 21¼")	N/A	N/A
	Overall width	580 mm (22¾")	630 mm (24¾")	690 mm (27¼")	690 mm (27¼")
	Length	640 mm (25¼")	720 mm (28¾")	840 mm (33")	870 mm (34¾")
	Length with starter kit	N/A	660 mm (26")	N/A	N/A
	Turning diameter	1440 mm (56¾")	1600 mm (63")	1900 mm (75")	2000 mm (78¾")
	Maximum user height	1050 mm (41¼")	1150 mm (45¼")	1500 mm (59")	1800 mm (70¾")

Size 0, 1, 2 & 3					
		mm (inch), kg (lb)			
		Size 0	Size 1	Size 2	Size 3
					
	Brake grip distance	74 mm (3")	74 mm (3")	74 mm (3")	74 mm (3")
	Static stability downhill	15°	15°	15°	15°
	Static stability uphill	7°	7°	7°	7°
	Static stability sideways	4.5°	4.5°	4.5°	4.5°
	Folded length	650 mm (25¾")	720 mm (28½")	860 mm (34")	1000 mm (39½")
	Folded width	580 mm (22¾")	620 mm (24½")	690 mm (27¼")	690 mm (27¼")
	Folded height	300 mm (12")	300 mm (11¾")	270 mm (10¾")	290 mm (11½")

Size 0, 1, 2 & 3

	mm (inch), kg (lb)			
	Size 0	Size 1	Size 2	Size 3
 Mass device	6.1 kg (13.4 lb)	6.4 kg (14.1 lb)	6.8 kg (15 lb)	11.3 kg (25 lb)
 Maximum user mass	30 kg (66 lb)	30 kg (66 lb)	45 kg (99 lb)	80 kg (176 lb)
 Maximum rated load	38.1 kg (84 lb)	38.4 kg (84.7 lb)	53.8 kg (118.6 lb)	93.3 kg (205.7 lb)
 Small castor (width /diameter)	W=28 mm/ D=140 mm (W=1"/ D=5½")	W=28 mm/ D=140 mm (W=1"/ D= 5½")	W=28 mm/ D=140 mm (W=1"/ D= 5½")	W=28 mm/ D=140 mm (W=1"/ D=5½")
 Large castor (width /diameter)	W=30 mm/ D=200 mm (W=1¼"/ D=8")	W=30 mm/ D=200 mm (W=1¼"/ D=8")	W=30 mm/ D=200 mm (W=1¼"/ D=8")	W=30 mm/ D=200 mm (W=1¼"/ D=8")
 Handgrip width	25 mm (1")	32 mm (1¼")	32 mm (1¼")	32 mm (1¼")
 Handgrip length	190 mm (7½")	190 mm (7½")	190 mm (7½")	250 mm (9¾")

Size 1 AT, 2 AT & 3 AT

Size 1 AT



Size 2 AT



Size 3 AT



Size 1 AT



<1190 mm
<46³/₄"



<30 kg
<66 lb

Size 2 AT



<1540 mm
<60³/₄"



<45 kg
<99 lb

Size 3 AT




<1840 mm
<72¹/₂"




<80 kg
<176 lb


Size 0, 1, 2 & 3




Size 0



Size 1



Size 2



Size 3

