

Bariatric Sentra Extra-Heavy-Duty Wheelchair

HCPCS: K0007,K0007



FIGURE A

- Weight Capacity: 500 lbs
- Carbon Steel frame with silver vein finish is attractive, and easy-to-maintain
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Reinforced side frame, caster journals and dual cross bars
- Dual-position front fork
- Dual-position front frame
- Composite Mag-style wheels are lightweight and maintenance free
- Hemi-adjustable caster is optional (Part #STDS1074)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Padded armrests provide added user comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD20DDA-ELR	SentraHD,20"DetDeskArmElvLeg,1/cs
STD20DDA-SF	SentraHD,20"DetDeskArm,Swg.Ft,1/cs
STD22DDA-ELR	SentraHD,22"DetDeskArmElvLeg,1/cs
STD22DDA-SF	SentraHD,22"DetDeskArm,Swg.Ft,1/cs
STD24DDA-ELR	SentraHD,24"DetDeskArmElvLeg,1/cs
STD24DDA-SF	SentraHD,24"DetDeskArm,Swg.Ft,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	54 lb 56 lb 66.25 lb
Product Dimensions	42" (L) w/ riggings x 29" 31" 32.5" (W) x 37" (H)
Seat Dimensions	20" 22" 24" (W) x 18" (D)
Weight Capacity	500 lb
Seat To Floor Height	19" - 21"
Carton Shipping Weight	60 lb 62 lb 72 lb
Folded Dimensions	42" (L) w/ riggings x 13" (W) x 37" (H)
Back Height	18"

ACCESSORIES	Description
WASL	Limb Support,Left,1/cs, CS

ACCESSORIES	Description
WASR	Limb Support, Right,1/cs, CS
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS802	Anti-Tipprs,w/Wheels,Adj,1pr/bx, BX